



Medium Duration Isometrics A



STRENGTH - GPP PROTOCOL

*Click workout title for video demonstration

*Perform entire workout as a circuit for 2-3 rounds

*Each movement is to be performed for the entire duration without any breaks

*Maintain nasal breathing during entire isometric

No.	Exercise	Load	Volume	Rest
A.	BB Righth Leg Split Squat	30-50%	30-60 seconds	1-2 min.
B.	DB Bench Press	30-50%	30-60 seconds	1-2 min.
C.	DB Chest Supported Row	30-50%	30-60 seconds	1-2 min.
D.	BB Left Leg Split Squat	30-50%	30-60 seconds	1-2 min.
E.	DB Z-Press	30-50%	30-60 seconds	1-2 min.
F.	BB RDL	30-50%	30-60 seconds	1-2 min.
F.	DB Pullover	30-50%	30-60 seconds	1-2 min.

<https://www.masqtraining.com/>